

CONTEC

Stock Code:300869



LIGHT THE RHYTHM OF LIFE WITH TECHNOLOGY



Digital Video Laryngoscope CMS-GS1

Contec Medical Systems Co.,Ltd.

ADD:No.112 Qinhuang West Street, Economic & Technical Development Zone,
Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA

Tel: +86-335-8015430

Fax: +86-335-8015588

E-mail: cms@contecmed.com.cn

Website: <http://www.contecmed.com>

CONTEC reserves the right of product appearance, specifications and function modification and updating,
without prior notification.



Technical Specification

Five solution for the Intubation



Safety

- Ergonomic and durable design for long time usage
- High quality material with anti shatter protection for whole unit
- Disposable Blades avoid cross infection
- Anti- Fog Lens and Anti-Microbial Handle



Fast

- Fast intubation after short time training
- Fast install or remove disposable blade
- Lightweight and portable



High- Definition

- 3.5 inch anti- reflective full touch display with 640*480 resolution
- 2.0 Mega Pixel Camera provide 1280*720P photo and video
- Rotation angle for up-down $\geq 120^\circ$ and left-right $\geq 270^\circ$



Easy

- One shot button for Recording
- One cable connect to PC for data management
- Full touch screen with user friendly operation



Cost Effective

- Favorable cost for all kinds of disposable blades
- Rechargeable battery for at least 2 hours working time

Suitable for

- Difficult intubation of glottic opening, easily guide the medical staff to accurately perform endotracheal intubation for anesthesia or first aid.

Technical Specifications

Display	Monitor: LCD
	Resolution Ratio: 640X480RGB
	Aspect Ratio: 4:3
	Video refresh Rate: 30FPS
Data output	Easy to create and store files in JPEG or MP4 format. Connect to the computer via a USB micro cable.
Camera and image	Resolution Ratio 2.0M Pixels
	Field Angle: $\geq 60^\circ$
	Illuminance $\geq 2500\text{LUX}$
	Working distance: 20mm-100mm
	Color temperature $\geq 5000\text{K}$
Power and Battery	Type Rechargeable Lithium Battery
	Battery life Cycle: More Than 500 times
	AC100-240V, 50/60Hz, 50VA
	Charging times: less than 2hrs
Working Condition	Temperature: -10°C ~ $+50^\circ\text{C}$
Size and weight	Product Size: 260.51*100.9*107.47mm; Product Weight: 0.22kg
	Carton size: 340*250*195mm ; Weight: 1.93kg



Mult-standard power adapter

HD lens



G4 Adult large(6.0mm or more)



G3 Adult (6.0~9.0mm)



G2 Child (3.0~6.0mm)